## Reinigungshinweise

### Reinigung

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## Cleaning Instructions

### Cleaning

Successful experimental results can only be achieved by using clean apparatus. The following advice may be used for the massive and porous borosilicate glass 3.3:

Wash glassware as quickly as possible after use. During washing with water all parts of the article should be scrubbed with a soft brush. Avoid any abrasion of the glassware! After cleaning thoroughly rinse with distilled water. Drying of the soaked filters can be undertaken either in air or in a dry box at temperatures not exceeding 100° C.

If after mechanical cleaning for example filter-pores should still remain clogged, the thorough chemical cleaning is required. A few suggestions found generally useful follow:

Contamination	Solvent
Barium sulphate	Hot concentrated sulphuric acid (100° C) silver-chloride, hot ammonia solution or sodium hyposulfite
Copper / Iron Oxides	Hot hydrochloric acid plus potassium chlorate
Mercury residue	Hot nitric acid
Mercury sulphide	Hot aqua regia
Albumen	Hot hydrochloric acid or hot ammonia
Grease, oil	Carbon tetrachloride
Organic matters	Hot concentrated cleaning solution, or hot concentrated sulfuric acid plus a few drops of sodium or potassium nitrite
Glucose	Carefully heat to approx. 200° C with mixed acid (5% H <sub>2</sub> SO <sub>4</sub> and 1% HNO <sub>3</sub> )

Extensive rinsing with water must obviously follow.

Initial cleaning before first use can also improve the filtration performance. We therefore recommend mild cleaning with warm water before first use. The glass filters can be rinsed and dried as described above.

### Warning!

Hot, concentrated phosphoric acid and hot alkaline solutions attack the glass surface, they are unsuitable as cleaning agents. If they have to be filtered, an increase in poresize and reduced life of the fritted disc is unavoidable.