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## Safety Instructions for Borosilicate 3.3 Glassware

### Cleaning

Before the first use of glass-filters, they should be cleaned with warm, diluted (~10%) hydrochloric acid, followed by several rinses with distilled water to remove dust and loose particles. Glassfilters should be cleaned immediately after their use. The surfaces can be wiped with a soft squeegee or brush. Surface rinsing with a wash bottle or back-flushing under the water tap is often sufficient.

### Pressure and Vacuum

Even small scratches or internal abrasions of the glass surface can greatly reduce the vacuum- or pressure resistance. Filter pores may clog during filtration and cause a pressure drop. In any case 0,1 Mpa may not be exceeded. Never subject glassware to sudden pressure changes.

### Temperature Resistance

The maximum allowable service temperature is approx. 500 degrees Centigrade. To avoid permanent stress, the heating and cooling rate should not exceed 5 °C/min. Wet glass-filters should slowly be heated up to 80 degrees C and dried for one hour, before temperature is increased. Permanent stress can reduce the mechanical and thermal resistance. Never subject glassware to sudden temperature changes.

### Chemical Resistance

High concentrations of hydrofluoric acid or hot alkali solutions attack the glass with increasing temperature and can cause a deterioration in the filter and a destruction of the pore size.

### Safety Advice

When working with glassware always wear protective glasses and protective gloves to avoid injuries. Use a safety screen, -hood or similar protection when working with glassware subject to pressure or vacuum. Be prepared, that reagents may leak from a broken vessel. Please heed these guidelines in combination with the respective state-specific regulations for the use of glassware in the laboratory.